












Glove Comparison Chart

Glove material	Intended use	Advantages and disadvantages	Example Photos
Latex (natural rubber)	Incidental contact	<p>Good for biological and water-based materials. Poor for organic solvents. Little chemical protection. Hard to detect puncture holes. Can cause or trigger latex allergies</p>	
Nitrile	<p>Incidental contact (disposable exam glove) Extended contact (thicker reusable glove)</p>	<p>Excellent general use glove. Good for solvents, oils, greases, and some acids and bases. Clear indication of tears and breaks. Good alternative for those with latex allergies.</p>	
Butyl rubber	Extended contact	<p>Good for ketones and esters. Poor for gasoline and aliphatic, aromatic, and halogenated hydrocarbons.</p>	
Neoprene	Extended contact	<p>Good for acids, bases, alcohols, fuels, peroxides, hydrocarbons, and phenols. Poor for halogenated and aromatic hydrocarbons. Good for most hazardous chemicals.</p>	
Norfoil	Extended contact	<p>Good for most hazardous chemicals. Poor fit (Note: Dexterity can be partially regained by using a heavier weight Nitrile glove over the Norfoil/Silver Shield glove.)</p>	
Viton	Extended contact	<p>Good for chlorinated and aromatic solvents. Good resistance to cuts and abrasions. Poor for ketones. Expensive.</p>	
Polyvinyl chloride (PVC)	Specific use	<p>Good for acids, bases, oils, fats, peroxides, and amines. Good resistance to abrasions. Poor for most organic solvents.</p>	
Polyvinyl alcohol (PVA)	Specific use	<p>Good for aromatic and chlorinated solvents. Poor for water-based solutions.</p>	
<p>Stainless steel Kevlar Leather</p>	Specific use	<p>Cut-resistant gloves. Sleeves are also available to provide protection to wrists and forearms. (If potential for biological or chemical contamination: wear appropriate disposable gloves on top of your cut-resistant gloves and discard after use).</p>	  
Cryogenic Resistant Material Leather	Specific use	<p>For use with cryogenic materials. Designed to prevent frostbite. Note: Never dip gloves directly into liquid nitrogen.</p>	