Level I

FUME HOOD SHUT-DOWN PROCEDURES

University of British Columbia Plant Operations

DO NOT USE THIS FUME HOOD

FOR REGIII.	AR MAINTENANCE WORK WHERE THE EXHAUST SYSTEM IS
TOR REGUE	STACKED ON THE ROOF AND SYSTEM IS INTACT
Planned Shut-	-Down Date: Time: AM/PM
DONE (√)	INSTRUCTIONS
` '	Ensure all containers are tightly capped.
	All <i>extremely hazardous</i> materials, such as compressed gas cylinders,
	must be removed from hood.
	Shut down all heat sources.
	Contact the Radiation Protection Office at 822-7072. Monitor the hood
	for radioactive decontamination where radioisotopes have been used and decontaminate if necessary.
	Ensure that no other lab equipment or apparatus is vented into the hood during shutdown.
	Lower the sash and attach this completed form to the front face.
	Advise all personnel in the area of the planned work.
Environment at 822-2029 to develop safe work procedures. I have prepared this good for REPAIR/MAINTENANCE as per the instructions given above. T be signed by the Hood User	
SIGNATURE	OF FUME HOOD USER: DATE:
CAUTION!!!	
	DO NOT REMOVE THIS TAG
	OR USE THIS FUME HOOD
	UNLESS SIGNED OFF BY THE TRADESPERSON
I have completed the REPAIR/MAINTENANCE and this fume hood may now be used. To be signed by the tradesperson.	
SIGNATURE	OF TRADESPERSON:DATE:
Blue	For additional copies, please contact HS&E 822-2029